

## DUFTDOS-DS

### FUNCTION DESCRIPTION

Peristaltic dosing pump for the dosing of fragrances into the steam pipe of steam rooms with intensity controller.

The dosing performance is adjusted by the pause time in between the injections by a rotary knob at the housing. The dosing time is set inside the housing. Dosing performance max. 3 l/h. Optionally the pump can be controlled by a push button from the steam room or by a non-volt contact.

### SCOPE OF DELIVERY

- Controller DS with rotary knob, installed in housing
- Peristaltic pump SA app. 3 l/h
- Dosing hose 3,2 x 1,6 Viton with hose holder
- Roller SA yellow
- Power supply cable with Euro plug
- Suction hose PE 4 x 1 - 1,5 m

### TECHNICAL DATA

|                    |                   |
|--------------------|-------------------|
| Electrical supply  | 230 V-AC, 50 Hz   |
| Power supply       | 6 W               |
| Back pressure max. | app. 0,5 bar      |
| Dimensions (W/H/D) | 80 x 180 x 140 mm |
| Weight             | app. 0,6 kg       |

### ACCESSORIES / OPTIONS

- Console with bottle 1,5 l - Part no. 13892
- Console for 10 l canister - Part no. 14391
- Dosing set consisting of: dosing valve stainless steel  $\frac{3}{8}$ " - 4 x 1 and 3 m dosing tube PTFE 4 x 1 - Part no. 12596

