

## DUFTDOS-0-T

### FUNCTION DESCRIPTION

Peristaltic dosing pump for the dosing of fragrances into the steam pipe of steam rooms. For external control with 230 V-AC e.g. by the steam generator or by an external controller. The control signal must be interlocked with the steam production. With tests switch for priming and test function.

### SCOPE OF DELIVERY

- Peristaltic pump SA app. 3 l/h
- Dosing hose 3,2 x 1,6 Viton with hose holder
- Roller SA yellow
- Connection cable 3 x 0,75<sup>2</sup>
- Suction hose PE 4 x 1 - 1,5 m
- Test button

### TECHNICAL DATA

Electrical supply	220-230 V-AC, 50-60 Hz
Power supply	6 W
Back pressure max.	app. 0,5 bar
Dimensions (W/H/D)	80 x 180 x 140 mm
Weight	app. 0,6 kg

### ACCESSORIES / OPTIONS

- Console with bottle 1,5 l - Part no. 13892
- Console for 10 l canister - Part no. 14391
- Dosing set consisting of: dosing valve stainless steel  $\frac{3}{8}$ " - 4 x 1 and 3 m dosing tube PTFE 4 x 1 - Part no. 12596

